

EXAMINING THE EYE in VR

Claire Chen and Luci Feinberg



WHO ARE WE WORKING WITH?





Dr. Stanley Chang

Ophthalmologist at CUIMC/Edward S. Harkness Eye Institute

Visualization of Eye Conditions

Slit Lamp Microscope (SLM)



Visualization of Eye Conditions

Doctors use a stereoscope to render stereo (3D) images through binocular vision







Right eye

However



Stereoscopes are not convenient to use

- Currently, observers and students can only view the 2D images through a single ocular.
- Not desirable for remote education and telemedicine.
- Doctors can't scale or move images freely for closer examination.
- Doctors can't compare and contrast different images of the same eye over time.

Visualization of Eye Conditions in VR



Goal 1: Render stereo (3D) images through binocular vision in Quest 3







Right eye

Visualization of Eye Conditions in VR

Test: creating different layers for each eye and displaying differently colored images for each.

Left eye layer





Right eye layer

Visualization of Eye Conditions in VR




Milestone 1: Rendered stereo (3D) still images in Quest 3, enabling scaling and ray grabbing



Feedback


- In addition to the images, display patient information and timestamp
- A library where we can sort and select the images

Enable the same features for stereo videos




DEMO






TECHNICAL IMPLEMENTATION



PLANS FOR THE UI




Viewer pane, only able to view 3 at a time

Centralized library





Viewer pane, view multiple images at a time

Centralized library



NEXT STEPS





THANK YOU

QUESTIONS?